# PREempt®

## WIPES

For disinfection of hard non-porous surfaces in healthcare establishments such as compound pharmacy labs, laboratories, clean rooms, spas, salons, and the aviation and mold remediation industries.

## Reference Sheet

### EFFICACY CLAIMS

Germicidal activity of this product was determined in accordance with the Canadian General Standards Board's standard CAN/CGSB-2.161-97.

#### TUBERCULOCIDAL: 3 MINUTES

*Mycobacterium terrae* (ATCC 15755)

#### VIRUCIDAL: 3 MINUTES

Proven effectiveness against the Poliovirus Type 1, Sabin strain type 1 (ATCC VR-192) which allows for a Broad-Spectrum Virucide claim against most enveloped and non-enveloped viruses.

HIV-1 Human Immunodeficiency virus (HIV), Strain HTLV-IIIB (HIV-1)

Feline Calicivirus, F9 Strain (ATCC VR-782), as a surrogate for Norwalk and Norwalk-like viruses

Rhinovirus serotype 14 (ATCC VR-284)

Canine Parvovirus (CPV), the Cornell strain (ATCC VR-2017)

Influenza Virus PR8 strain (ATCC VR-95)

Rotavirus Wa strain (ATCC VR-2018)

Murine Norovirus type 1 (MNV-1) strain S99

Human Coronavirus 229E (ATCC VR-740)

Human Parainfluenza Virus 3 (ATCC VR-92)

Adenovirus type 5 (ATCC VR-5)

This product has demonstrated effectiveness against Poliovirus and is expected to inactivate all Influenza A viruses including 2009 (H1N1) pandemic Influenza A virus

#### BACTERICIDAL: 3 MINUTES

*Pseudomonas aeruginosa* (ATCC 15442)

*Staphylococcus aureus* (ATCC 6538)

*Salmonella enterica* (ATCC 10708)

*Methicillin-resistant Staphylococcus aureus (MRSA)* (ATCC 29247)

*Vancomycin-resistant Enterococcus faecalis (VRE)* (ATCC 51299)

*Escherichia coli O157:H7* (ATCC 43888)

*Enterobacter aerogenes* (ATCC BAA-2356)

#### FUNGICIDAL: 3 MINUTES

*Trichophyton mentagrophytes* (ATCC 9533)

*Candida albicans* (ATCC 10231)

#### BROAD-SPECTRUM SANITIZING: 30 SECONDS

*Pseudomonas aeruginosa* (ATCC 15442)

*Staphylococcus aureus* (ATCC 6538)

*Salmonella enterica* (ATCC 10708)

*Methicillin-resistant Staphylococcus aureus (MRSA)* (ATCC 29247)

*Vancomycin-resistant Enterococcus faecalis (VRE)* (ATCC 51299)

*Klebsiella pneumoniae* (ATCC 13882)

*Acinetobacter baumannii* (ATCC 19606)

*Escherichia coli O157:H7* (ATCC 43888)

*Campylobacter jejuni* (ATCC 33560)

### DISINFECTANT USE DIRECTIONS

Heavily soiled surfaces require cleaning prior to disinfection.

**TO DISINFECT NON-CRITICAL, PRE-CLEANED INSTRUMENTS AND MEDICAL DEVICES:** Instruments must be thoroughly pre-cleaned to remove excess organic debris, rinsed, and dried. Wipe the instrument with PREempt Wipes and allow to remain wet for 3 minutes. Wipe with a clean, damp cloth or paper towel and allow to air dry.

**DISINFECTION OF NON-POROUS HARD SURFACES** such as the exterior of fume hoods, tabletops and other environmental surfaces, etc. Apply to surface with disposable wipe. Ensure surface remains wet for 3 minutes.

#### TO PULL WIPES THROUGH LID, SIMPLY...

1. Remove lid.
2. Remove liner from the canister.
3. Pull up corner of first wipe from centre of roll and push 1 to 2 cm, ½" to ¾" through "X" in lid.
4. Replace lid.
5. Lift lid flip top and pull up wipes one at a time.

To close between uses to prevent drying, simply keep lid flip top down and ensure wipe corner is inside the raised circle in lid.

**BROAD-SPECTRUM SANITIZING ON ENVIRONMENTAL SURFACES (>99.9%):** Apply to surface with disposable wipe, allow to remain wet for 30 seconds. Wipe dry. No rinse required. Food contact surfaces require rinsing with potable water.

#### SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION

**AGAINST HIV** (Human Immunodeficiency Virus) on objects and surfaces soiled with blood/body fluids. This product is intended for use against HIV only in those settings where the virus would be expected to be encountered, such as settings where contamination by blood or body fluids is likely.

#### CLEANING AND DISINFECTING SURFACES OF BLOOD AND BODY FLUIDS:

Gloves should be worn. Remove excess blood and fluid with absorbent materials. **Clean contaminated area:** Wipe the surface with disposable PREempt Wipes. Ensure all blood/body fluids are thoroughly cleaned from surface/objects before starting disinfection. **Disinfect contaminated area:** Apply to surface, allow surface to remain wet for 3 minutes.

**PERSONAL PROTECTION:** Disposable gloves, gowns, facemasks, or eye coverings as appropriate, must be worn during all cleaning of body fluids, blood and decontamination procedures.

**DISPOSAL OF INFECTIOUS MATERIAL:** Products contaminated with blood or body fluids should be disposed of according to Federal, Provincial, and local regulations for infectious waste disposal.

### PRECAUTIONARY STATEMENTS

No Personal Protective Equipment is required when handling this product. **KEEP OUT OF REACH OF CHILDREN.** Avoid storage at elevated temperatures. Avoid contamination of food in the application and storage of this product. Do not mix with other cleaning or disinfecting products.

**FIRST AID:** If in contact with eyes, flush immediately and thoroughly with water. Call a physician if irritation persists. If ingested, do not induce vomiting. Call a physician or a poison control center immediately. TAKE CONTAINER, LABEL OR PRODUCT NAME AND DIN WITH YOU WHEN SEEKING MEDICAL ATTENTION.

**DISPOSAL:** For information on the disposal of unused, unwanted product and the clean up of spills, contact the Provincial Regulatory Agency or the Manufacturer. Rinse the container with water. Dispose of container in accordance with municipal, provincial and federal regulations. Do not reuse containers.

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