

What's the Difference between the AS SUPER 4 Digital and the ES-130 / IC-1107?

Features	AS SUPER 4 Digital	ES-130 / IC-1107
3 channel		X
4 channel	X	
Digital display	X	
Rotary Dial (0-10)		X
Digital Increments (+/- .1 mA up to 20mA max/channel)	YES	NO
Includes battery	AAA	9-volt
Adaptor Included?	YES	NO
TENS	NO	YES
Programmable	Yes - 37 Programs included	NO
Nogier Frequencies	YES	NO
Free Han Stimulation	YES	NO
Frequency Modulation (FM)	YES	NO
Burst Programs	YES	NO
TSEA Programs	YES	NO
Frequency Range	1-120Hz	1-500Hz
Warranty (Device only, does not include accessories)	2 YEAR	1 YEAR