

Introducing the next generation of a clinical icon.



Open bell versatility.

The pediatric side converts to a traditional open bell by replacing the single-piece diaphragm with a non-chill rim.

Introducing the 3M™ Littmann® Classic III™ Stethoscope

Right out of the box, it lives up to its name.

The Littmann Classic III stethoscope is the latest version of the stethoscope that helps millions of medical professionals achieve their best. It brings new design, materials, and technology to the Littmann stethoscope series trusted by clinicians around the world. The new Classic III stethoscope offers you:



High acoustic sensitivity for performing general physical assessments.



Tunable diaphragms on both the adult and pediatric sides of the chestpiece.



Single-piece diaphragm that's easy to attach and easy to clean because the surface is smooth without crevices.



Next-generation tubing with improved resistance to skin oils and alcohol for longer tubing life. Less likely to pick up stains.

Using a Classic III stethoscope demonstrates what's important to you as a medical professional: quality, achievement, and providing the best possible care. Every patient. Every day.

Here's what makes the

3M™ Littmann® Classic III™ Stethoscope stand out.

Features

Tunable technology is a 3M invention that lets you shift between hearing high and low frequency sounds by simply adjusting how much you press on the chestpiece.

Soft-sealing eartips give you a comfortable fit and an excellent acoustic seal; they snap tightly onto the eartubes for safety.

Open bell stays clear of dirt and debris by keeping it covered with the small diaphragm.

Adjustable headset lets you easily set the tension for head size and comfort just by squeezing or pulling apart the eartubes.

No natural rubber latex or phthalate plasticizers are used in the tubing or any other component—to help protect human health and the environment.

Made in the USA for consistent quality. We apply rigorous process and quality controls to every Littmann stethoscope. You'll always get a durable, dependable, clinical tool that will serve you for years to come.

Classic III stethoscope specifications

Length	69cm / 27.2 inches
---------------	--------------------

Weight	150g
---------------	------

Acoustic performance	7
-----------------------------	---

Warranty	5 years
-----------------	---------

CHESTPIECE

Construction	Two-sided, adult and pediatric
---------------------	--------------------------------

Material	Stainless steel
-----------------	-----------------

Diaphragm type	Tunable, single-piece
-----------------------	-----------------------

Diaphragm diameter	Adult: 4.3cm / 1.7 in Pediatric: 3.3cm / 1.3 in
---------------------------	--

TUBING

Construction	Single lumen (single acoustic channel within tube)
---------------------	---

Material	Not made with natural rubber latex or phthalate plasticizers
-----------------	---

HEADSET

Material	Wide diameter aerospace alloy / anodized aluminum
-----------------	--