

## Jaybird Comparison Chart

	Fabric	Rigid or Stretch	Bleached	Latex Free
<b>Jaybird Poly Cotton Blend Tape</b>	40% Polyester 60% Cotton blend	Rigid	No	Yes
<b>Jaybird One Tape</b>	50% Polyester 50% Cotton	Rigid	Yes	Yes
<b>Jaybird Pro-White Tape</b>	40% Polyester 60% Cotton	Rigid	Yes	Yes
<b>Jaybird 2550 Lightweight Stretch Tape</b>	20% Polyester, 75% Cotton & 5% Spandex	Stretch	No	Yes
<b>Jaybird Victor Rigid Strapping Tape</b>	100% Rayon, dyed textile	Very rigid (equivalent to Leukotape P)	No	No

	Computerized Unwind Tension	Serrated Edge Cut	Hand Tearable	Tensile Str (lbs)
<b>Jaybird Poly Cotton Blend Tape</b>	Yes	Yes	Yes	45
<b>Jaybird One Tape</b>	Yes	Yes	Yes (bidirectional)	50
<b>Jaybird Pro-White Tape</b>	Yes	Yes	Yes (bidirectional)	43
<b>Jaybird 2550 Lightweight Stretch Tape</b>	No	No	Yes	31
<b>Jaybird Victor Rigid Strapping Tape</b>	No	Yes	Yes	72