

The logo for OrthoCanada, featuring the word "ORTHOCANADA" in white, uppercase, sans-serif font on a red rectangular background.

The Martindale rating and what it means to you:

SEERS Medical builds its tables using components of the highest quality, that is reflected in the Martindale rating for vinyl upholstery. The Martindale rating is the measuring scale for upholstery products and requires the use of calibrated tools to assess the durability of materials. Simply put, the higher the rating, the more durable the fabric. Industry norm is at a Martindale rating of 30 000 rubs, the Mundial vinyl selected by SEERS Medical is rated at 100 000, indicating that it is a highly durable fabric that will last a long time.

[Polartec Fabric Testing - Martindale Abrasion Test \[YouTube Video\]](#)

The logo for OrthoCanada, featuring the word "ORTHOCANADA" in white, uppercase, sans-serif font on a red rectangular background.

For more information, please see our contact information below:

1-800-561-0310 | info@orthocanada.com | www.orthocanada.com