

# NEEDLE COMPARISON



NEEDLE COMPARISON	JAPANESE	KOREAN	CHINESE
	SEIRIN	DONG BANG	HWA-TO
SILICONE COATED SHAFT	✓ (JUST THE TIP)	✓	
TABS	ULTRASONIC BINDING	✓	✓
SMOOTH POLISH FINISH	✓		
METAL HANDLE	PLASTIC HANDLE	✓	✓
BULK PACKAGING	LIMITED SIZES	✓	
FULL RANGE OF SIZES	WIDEST RANGE OF GAUGE & LENGTHS	LONGEST IS 60MM	LIMITED GAUGE SIZES
PRESS NEEDLES	✓	✓	ONLY IN CARBO