

Specifications

	Pro TENS	Pro TENS+	Pro TENS / EMS
Display	Analog	LCD	LCD
Channels	Dual, isolated between channels	Dual, isolated between channels	Dual, isolated between channels
Waveform	Asymmetrical bi-phasic square pulse	Asymmetrical bi-phasic square pulse	Asymmetrical bi-phasic square pulse
Amplitude	0-80 mA	0-80 mA	0-80 mA
Pulse Rate	2Hz-120Hz	1Hz-160Hz adjustable, 1Hz/step (1-20Hz), 5Hz/step (20-160Hz)	1Hz-160Hz adjustable, 1Hz/step (1-20Hz), 5Hz/step (20-160Hz)
Pulse Width	40µs-260µs	50µs-260µs adjustable, 10µs/step	TENS: 50µs-260µs adjustable, 10µs/step EMS: 250µs fixed
Modes	B, C & M	B, C & M, SDR & SDW	TENS: B, C & M SDR & SDW EMS: S & A
Preset Programs	No	No	No
Program Lock	No	Yes	Yes
Compliance Meter	No	Yes	Yes
Timer	30, 60 min & C	15, 30, 60 min & C	15, 30, 60 min & C
Battery	9V	9V	9V
Low Battery Indicator	Yes	Yes	Yes
Size (L x W x D)	3.7" x 2.4" x 1"	4.3" x 2.4" x 0.98"	4.3" x 2.4" x 0.98"
Weight	120g with battery	140g with battery	140g with battery
Warranty	1 year	1 year	1 year
MDALL	108683	108683	108683
Code	OC-32035	OC-32036	OC-32037