

SEERS Hi-Lo Adult Size Change Table

SKU V-220-B | V-230-B

Upholstery

Antimicrobial Additive (Sanitized®)

Seamless Upholstery with Easy-Clean Backing

22oz PVC, Martindale 150,000 Rubs, CRIB 5 Fire Resistance

Frame

Antimicrobial Powder Coating

3mm Thick Steel Tubing

Robotically Welded Interlocking Joints

Self-lubricating Precision Bearings

4" Elevation Offset

18" - 39" Electric Height Adjustment

440lbs Lifting Capacity, 1,113lbs Static Capacity

Optional Side Rails and Side Rail Cushions


Exceptional Warranty

Lifetime - Frame

5 years - Electronics, Mechanics, Components




V-220-B model shown with optional side rails and side rail cushions


SCAN FOR
MORE DETAILS



V-220-B



V-230-B



440
Lifting
Capacity (lbs.)

Accessible

Anti-
Microbial

Freight
Required

Meets Canadian electrical & load safety standards

CSA®

PS
Certified

Change Table

Change Table with Backrest

Dimensions	27.5"W x 72"L x 18"-39"H	27.5"W x 76.8"L x 18"-39"H
Backrest Length	N/A	27.6"
Backrest Tilt	N/A	-20°, +75°
Foot Section Length	N/A	49.2"
1-Hand Release	N/A	Yes
Foam	1.6" medium density	
Upholstery	Waterproof, antimicrobial & fire retardant Mundial vinyl	
Martindale Abrasion Rating	150 000 rubs	
Lifting Capacity	440lbs	
Travel Time	45 sec	
Castor System	Individually locking castors	
Side Rails & Cushions	Optional	
Clearance Beneath Frame	7 1/8"	
Paper Roll Holder	Optional	
Safety Marking	QPS/CSA	
Warranty	Lifetime frame, 5 yrs components	
SKU	V-220-B	V-230-B