



PRODUCT DATA SHEET

#2550LF Jaybird® LST LIGHTWEIGHT ADHESIVE STRETCH TAPE

DESCRIPTION

- ✦ Off White Cotton/Polyester/Spandex® blended fabric tape
- ✦ #2550 is coated with an aggressive zinc oxide hypoallergenic white synthetic rubber adhesive formulation
- ✦ Non-sterile, for external usage and disposable.
- ✦ Jaybird branded two (2) color core w/parchment paper outer wrap.
- ✦ Hand Tearable
- ✦ 7½ yard length roll stretched
- ✦ Long stretch product with memory
- ✦ Stable Adhesive Coating when exposed to heat
- ✦ Store Cool & Dry for maximum performance
- ✦ Made in USA by Jaybird & Mais, Inc.
- ✦ Not made with Natural Rubber Latex
- ✦ 100% Latex Free

PHYSICAL PROPERTIES*

TEST/REFERENCE	AVERAGE VALUE
Tensile Strength (PSTC-31)	31 lbs/linear inch
Adhesion 180° Peel to SS (PSTC-101)	20-ounces/linear inch

* Pressure Sensitive Tape Council methods used when applicable. Testing represents typical performance values (averages).

STANDARD SIZE & PACKAGING

PRODUCT NO	SIZE	SIZE (METRIC)	ROLLS/CASE
2550-1075LF	1" x 7½ yards	2,5cm x 6,9m	48
2550-1575LF	1½" x 7½ yards	3,8cm x 6,9m	32
2550-2075LF	2" x 7½ yards	5,0cm x 6,9m	24
2550-3075LF	3" x 7½ yards	7,5cm x 6,9m	16

RECOMMENDED STORAGE

The above product should be stored in a cool dry environment for best results of product performance. Average recommended temperature for storage: 60° to 72° Fahrenheit (19.4° to 22.8° Centigrade)

Technical information provided by Jaybird & Mais, Inc. is based on experience and/or tests believed to be reliable, but their accuracy is not guaranteed and the results may not be relevant to every user's application. For this reason, Jaybird & Mais, Inc. does not accept responsibility or liability, direct or consequential, arising from reliance upon any information provided and the user should determine the suitability of the products for their intended use.

