



PRODUCT DATA SHEET

#40 Polyester Cotton Blend RIGID WHITE ADHESIVE TAPE

DESCRIPTION

- ✦ Unbleached cotton and polyester woven fabric based adhesive tape
- ✦ Proprietary White 100% Synthetic Rubber Adhesive Formulation
- ✦ Made in USA by Jaybird & Mais, Inc.
- ✦ Computerized Unwind Tension
- ✦ Serrated Edge Cut (SEC)
- ✦ Non-sterile, for external usage and disposable
- ✦ Hand Tearable
- ✦ Stable Adhesive Coating when exposed to heat
- ✦ Store Cool & Dry for maximum performance
- ✦ 100% Latex Free
- ✦ Not Made with Natural Rubber Latex

PHYSICAL PROPERTIES*

| Test/Reference | Average Value |
|-------------------------------------|-----------------------|
| Tensile Strength (PSTC-31) | 45 lbs. @ 1½ width |
| Adhesion 180° Peel to SS (PSTC-101) | 29-ounces/linear inch |
| Adhesion 180° Peel to Backer (101) | 24-ounces/linear inch |
| Thickness | 12 mils |

STANDARD SIZE & PACKAGING

| Product No | Size | Size (Metric) | Rolls/Cs |
|------------|----------------|---------------|----------|
| 40-1015WH | 1" x 15 yards | 2,5cm x 13,7m | 48 |
| 40-1515WH | 1½" x 15 yards | 3,8cm x 13,7m | 32 |
| 40-2015WH | 2" x 15 yards | 5,0cm x 13,7m | 24 |
| 40-1010WH | 1" x 10 yards | 2,5cm x 13,7m | 200 |

RECOMMENDED STORAGE

The above product should be stored in a cool dry environment for best results of product performance. Average recommended temperature for storage: 60° to 72° Fahrenheit (19.4° to 22.8° Centigrade)

Technical information provided by Jaybird & Mais, Inc. is based on experience and/or tests believed to be reliable, but their accuracy is not guaranteed and the results may not be relevant to every user's application. For this reason, Jaybird & Mais, Inc. does not accept responsibility or liability, direct or consequential, arising from reliance upon any information provided and the user should determine the suitability of the products for their intended use.

