

PRODUCT DATA SHEET

50 Series

FOAM UNDERWRAP TAPE (NON-ADHESIVE)



DESCRIPTION

- #50 is Polyurethane foam used as a pre-tape base that is applied over the skin and under our adhesive tape to protect the skin from prolonged adhesive contact.
- Highly porous breathable taping base
- This product is designed for topical (external human) purposes only.
- Standard Color (Tan) 50-27530
- Product does not contain adhesive
- Hand tearable
- Plain Core & Plain White Box
- Made in USA
- Does not contain Natural Rubber Latex

PHYSICAL PROPERTIES*

TEST/REFERENCE	AVERAGE VALUE
Tensile Strength (PSTC-31)	N/A
Adhesion 180° Peel to SS (PSTC-101)	0-ounces/linear inch
Ignition Temperature	800° - 850°
Extinguisher Media	Carbon Dioxide, Dry Powder
Thickness	.025" (+0.013" +/- 0.003") (0.635mm)

* Pressure Sensitive Tape Council methods used when applicable. Testing represents typical performance values (averages)

STANDARD SIZE & PACKAGING

Product No	Size	Rls/Cs	Product No	Size	Rls/Cs
50-27530	2¾" x 30 yards (7,0cm x 27,4m)	48	50-BLACK	2¾" x 30 yards (7,0cm x 27,4m)	48
50-RED	2¾" x 30 yards (7,0cm x 27,4m)	48	50-GOLD	2¾" x 30 yards (7,0cm x 27,4m)	48
50-BLUE	2¾" x 30 yards (7,0cm x 27,4m)	48	50-GREEN	2¾" x 30 yards (7,0cm x 27,4m)	48
50-ORANGE	2¾" x 30 yards (7,0cm x 27,4m)	48	50-PINK	2¾" x 30 yards (7,0cm x 27,4m)	48
50-BRGREEN	2¾" x 30 yards (7,0cm x 27,4m)	48			

RECOMMENDED STORAGE

The above product should be stored in a cool dry environment for best results of product performance. Average recommended temperature for storage: 60° to 72° Fahrenheit (19.4° to 22.8° Centigrade)

Technical information provided by Jaybird & Mais, Inc. is based on experience and/or tests believed to be reliable, but their accuracy is not guaranteed and the results may not be relevant to every user's application. For this reason, Jaybird & Mais, Inc. does not accept responsibility or liability, direct or consequential, arising from reliance upon any information provided and the user should determine the suitability of the products for their intended use.