



PRODUCT DATA SHEET

EX1 Jaybird® One PREMIUM RIGID WHITE ADHESIVE TAPE

DESCRIPTION

- ✦ Bleached white woven fabric based adhesive tape
- ✦ Increased thread count over EX25 Pro-White™
- ✦ Proprietary White 100% Synthetic Rubber Adhesive Formulation
- ✦ Made in USA by Jaybird & Mais, Inc.
- ✦ Computerized Unwind Tension
- ✦ Serrated Edge Cut (SEC)
- ✦ Non-sterile, for external usage and disposable
- ✦ Hand Tearable (bi-directional)
- ✦ Stable Adhesive Coating when exposed to heat
- ✦ Store Cool & Dry for maximum performance
- ✦ 100% Latex Free
- ✦ Not made with Natural Rubber Latex

PHYSICAL PROPERTIES*

TEST/REFERENCE	AVERAGE VALUE
Tensile Strength (PSTC-31)	50 lbs./linear inch @ 1½ width
Adhesion 180° Peel to SS (PSTC-101)	36-ounces/linear inch
Adhesion 180° Peel to Backer (PSTC-101)	29-ounces/linear inch
Thickness	11 mils

* Pressure Sensitive Tape Council methods used when applicable. Testing represents typical performance values (averages)

STANDARD SIZE & PACKAGING

Product No	Size	SIZE (METRIC)	ROLLS/CASE
EX1-1015W	1" x 15 yards	2,5cm x 13,7m	48
EX1-1515W	1½" x 15 yards	3,8cm x 13,7m	32
EX1-2015W	2" x 15 yards	5,0cm x 13,7m	24

RECOMMENDED STORAGE

The above product should be stored in a cool dry environment for best results of product performance. Average recommended temperature for storage: 60° to 72° Fahrenheit (19.4° to 22.8° Centigrade)

Technical information provided by Jaybird & Mais, Inc. is based on experience and/or tests believed to be reliable, but their accuracy is not guaranteed and the results may not be relevant to every user's application. For this reason, Jaybird & Mais, Inc. does not accept responsibility or liability, direct or consequential, arising from reliance upon any information provided and the user should determine the suitability of the products for their intended use.

