

VERSA-BATH™ HEATING UNIT

- **Versatile design** - use as a splinting bath and/or moist heat pack unit
- **Easy clean** seamless design with rounded corners
- Easy to read, adjustable **digital temperature control**
- **Energy efficient**, insulated, double wall design
- No mess removable lid doubles as drip tray
- Large capacity rack holds 4 standard and 2 oversize moist heat packs
- Concealed heating element will not corrode and **lasts much longer** than exposed elements



VERSATILE DESIGN

The Versa-Bath™ is the first multi-purpose water bath that can be used for conditioning Thermal Core Moist Heat Packs, or as a water bath for molding thermoplastic splinting and casting materials. The large capacity rack holds 4 standard and 2 oversize moist heat packs and can be easily removed for splinting/casting applications. The no mess, removable lid doubles as a drip tray to keep your cabinets and floors free of water when removing the packs or transporting splinting material.

ENERGY EFFICIENT

These units incorporate an insulated, double wall design and an easy to read, adjustable digital temperature control to ensure proper water temperature for both moist heat therapy and splinting applications. The Versa-Bath heats water quickly and maintains a consistent temperature resulting in lower energy consumption.

EASY TO CLEAN

The seamless design, with rounded corners along with the concealed heating element, makes cleaning a snap! Because the heating element is not exposed to water, it will not corrode, lasting much longer than exposed elements.

CSA CERTIFIED

All Versa-Bath Heating Units are built to ANSI Standards and certified by the Canadian Standards Association. The CSA mark on your heating unit ensures that it has been evaluated to meet or exceed the ANSI/UL/CSA consumer safety standards for use in the United States and Canada.

SPECIFICATIONS

- Versa-Bath 120V Unit
- 4 Std Thermal Core® Packs
- 2 Cervical Thermal Core® Packs

Product Weight: 14.25 lb

Materials: 304 Stainless Steel

Preset Temp Range: 140°-180° F
60°-82° C

Dimensions:

Internal: 14"W x 9"D x 10.5"H
35.5cm x 23cm x 27cm

External: 15"W x 10"D x 16"H
38cm x 25.5cm x 41cm